

FARNAX®-CODIP



Optimal mastitis-prevention after milking

- Strong disinfecting
- Maximum protection in difficult circumstances
- Faster sphincter closing



WHEN TO USE?

- In cases where there is a high risk of mastitis
- When there is a high infection pressure: humid climate, wet, dirty bedding,...

EFFICACY

Farnax®-CoDip gives an optimal protection against different types of mastitis, caused by bacteria:

- Environmental
- Clinical
- Subclinical

COMPOSITION


Farnax®-CoDip is a 2-compound teat dipping product, based on chlorine dioxide. The high amount of emollients (sorbitol and glycerol) results in a unique formula with maximum teat protection.

Thanks to its high covering properties, Farnax®-CoDip supports a faster sphincter closing and so reduces the time where bacteria could enter the teat.

INSTRUCTIONS FOR USE



- Mix both compounds and let react for 1 – 4 hours (after mixing the product is stable for at least 3 weeks)
- Apply the product with a dipping cup



Use biocides with precaution. Before any use, read the label and the information concerning the product. This product is not necessarily available or registered in every country.